

Product Datasheet

RPS6KA1 (Phospho-T359+S363) antibody (orb572475)

Description	Recombinant rabbit monoclonal antibody to RPS6KA1 (Phospho-T359+S363)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	A synthesized peptide derived from human Phospho-RPS6KA1 (T359+S363)
Target	RPS6KA1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q15418
Dilution Range	WB: 1:500-1:5000, IF: 1:20-1:200, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.