

Product Datasheet

FOXO3 (Phospho-S253) antibody (orb572491)

Description	Recombinant rabbit monoclonal antibody to FOXO3 (Phospho-S253)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A synthesized peptide derived from human Phospho-FOXO3 (S253)
Target	FOXO3
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	O43524
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.