

Product Datasheet

NFKBIA (Phospho-S32) antibody (orb572498)

Description	Recombinant rabbit monoclonal antibody to NFKBIA (Phospho-S32)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	A synthesized peptide derived from human Phospho-NFKBIA (S32)
Target	NFKBIA
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P25963
Dilution Range	WB: 1:500-1:5000, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.