

Product Datasheet

PRKCZ (Phospho-T560) antibody (orb572519)

Description	Recombinant rabbit monoclonal antibody to PRKCZ (Phospho-T560)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthesized peptide derived from human Phospho-PRKCZ (T560)
Target	PRKCZ
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q05513
Dilution Range	WB: 1:500-1:5000
Expiration Date	12 months from date of receipt.