

Product Datasheet

SMAD2 (Phospho-S255) antibody (orb572526)

Description	Recombinant rabbit monoclonal antibody to SMAD2 (Phospho-S255)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	A synthesized peptide derived from human Phospho-SMAD2 (S255)
Target	SMAD2
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q15796
Dilution Range	IF: 1:20-1:200
Expiration Date	12 months from date of receipt.