

Product Datasheet

CDC37 (Phospho-S13) antibody (orb572527)

Description	Recombinant rabbit monoclonal antibody to CDC37 (Phospho-S13)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthesized peptide derived from human Phospho-CDC37 (S13)
Target	CDC37
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q16543
Dilution Range	WB: 1:500-1:5000
Expiration Date	12 months from date of receipt.