

Product Datasheet

PRKDC (Phospho-S2056) antibody (orb572554)

Description	Recombinant rabbit monoclonal antibody to PRKDC (Phospho-S2056)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	A synthesized peptide derived from human Phospho-PRKDC (S2056)
Target	PRKDC
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P78527
Dilution Range	IF: 1:20-1:200
Expiration Date	12 months from date of receipt.