

Product Datasheet

RB1 (Phospho-S780) antibody (orb572555)

Description	Recombinant rabbit monoclonal antibody to RB1 (Phospho-S780)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP
Immunogen	A synthesized peptide derived from human Phospho-RB1 (S780)
Target	RB1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P06400
Dilution Range	IHC-P: 1:50-1:200, IF: 1:20-1:200, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.