

Product Datasheet

MAP3K3 antibody (orb572567)

Description	Recombinant rabbit monoclonal antibody to MAP3K3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	A synthesized peptide derived from human MAP3K3
Target	MAP3K3
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q99759
Expiration Date	12 months from date of receipt.