

Product Datasheet

CDC5L antibody (orb572569)

Description	Recombinant rabbit monoclonal antibody to CDC5L
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human CDC5L
Target	CDC5L
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q99459
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:50-1:200, IF: 1:20-1:200, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.