

Product Datasheet

HSP90AA1 antibody (orb572618)

Description	Recombinant rabbit monoclonal antibody to HSP90AA1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IP, WB
Immunogen	A synthesized peptide derived from human HSP90AA1/HSP90AB1
Target	HSP90AA1/HSP90AB1
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P07900, P08238
Dilution Range	WB: 1:500-1:5000, IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.