

Product Datasheet

HCV NS4 (1691-1710)(1712-1733)(1921-1940) mosaic 34 kDa (orb572689)

Description	The protein contains immunodominant regions of HCV NS4, amino acids (1691-1710)(1921-1940)(1712-1733...
Form/Appearance	liquid (Supplied in Tris-HCl 50 mM pH 8.0, Sodium Chloride 60 mM, Glutathione 10 mM, Glycerol 50%)
Storage	store at -20 °C
Note	For research use only
Application notes	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Clonality	Recombinant
MW	Theoretical MW: about 34 kDa
Expiration Date	6 months from date of receipt.