

Product Datasheet

HIV-1 M gp41 (537-811) 58 kDa (orb572703)

Description	The protein contains immunodominant region of HIV-1 group M gp41, amino acids 537-811. The protein i...
Form/Appearance	liquid (Supplied in Imidazole 20 mM pH 7.0, Sodium Chloride 150 mM, Urea 8M, DTT 10 mM)
Storage	store at 4 °C
Note	For research use only
Application notes	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Clonality	Recombinant
MW	Theoretical MW: about 58 kDa
Expiration Date	6 months from date of receipt.