



Product Datasheet

DLinDMA (orb572878)

Description	DLinDMA is a key lipid component of stable nucleic acid lipid particles as a benchmark.
Storage	2 years in -20°C Powder, 2 weeks in 4°C in DMSO, 6 months in -80°C in DMSO
Note	For research use only
Formula	$C_{41}H_{77}NO_2$
MW	616.06
CAS Number	871258-12-7
SMILES	<chem>CCCCC/C=C/C/C=C\CCCCCCCCOCC(OCCCCCCCC/C=C/C/C=C\CCCC)C N(C)C</chem>
Expiration Date	12 months from date of receipt.



Chemical
Structure
(871258-12-7)



CAS: 871258-
12-7