

Product Datasheet

Super-TDU (orb573388)

Description	Super-TDU is an inhibitory peptide targeting YAP-TEADs interaction.
Note	For research use only
Formula	C ₂₃₇ H ₃₇₀ N ₆₆ O ₆₉ S
MW	5279.94
CAS Number	[1599441-71-0]
SMILES	<chem>CC[C@@H](C)[C@@H](C(=O)NCC(=O)NCC(=O)N[C@@H](CO)C(=O)NC(=O)N[C@@H](CC(=O)N)C(=O)N1CCC[C@H]1C(=O)N[C@@H](CCCCN)C(=O)N[C@@H]([C@H](C)O)C(=O)N[C@H](C)C(=O)N[C@@H](CC(=O)N)C(=O)N[C@@H](C(C)C)C(=O)N2CCC[C@H]2C(=O)C(=O)[C@H](CCC(=O)N)NN[C@@H]([C@H](C)O)C(=O)N[</chem>
Expiration Date	12 months from date of receipt.