



Product Datasheet

AEBP1 antibody (orb573858)


Description	Rabbit polyclonal antibody to AEBP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human AEBP1
Target	AEBP1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	Batch dependent within range: 100 µl at 0.5 - 1 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VPPEKTKDKGKKGGKDKGPKVPKESLEGSPRPPKKGKEKPPKATKKPKKEK
Clonality	Polyclonal
MW	131 kDa
Uniprot ID	Q14113
NCBI	NP_001120
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of HepG2 cells lys...



Immunohistochemical staining of human in...