

Product Datasheet

OR13C9 antibody (orb573949)

Description	Rabbit polyclonal antibody to OR13C9
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human OR13C9
Target	OR13C9
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8NGS9
Expiration Date	12 months from date of receipt.