


Product Datasheet

ATP1A1 (phospho-Tyr10) antibody (orb5746)

Description	ATP1A1 (phospho-Tyr10) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Na,K-ATPase alpha-1 around the phosphorylation site of tyrosine 10 (DK(p-Y)EP)
Target	Phospho-ATP1A1 (Tyr10)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	P05023
Dilution Range	ICC=1:100, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human bl...