

Product Datasheet

TBC1D10B antibody (orb575132)

Description	Rabbit polyclonal antibody to TBC1D10B
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rat, Yeast, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human TBC1D10B
Target	TBC1D10B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QKQEKEREKQEKEREKQEKERQKQEKKAQGRKLSLRRKADGPPGPHDGGD
Clonality	Polyclonal
MW	61kDa
Uniprot ID	Q9H8Z2
NCBI	NP_056342
Expiration Date	12 months from date of receipt.