

Product Datasheet

Estrogen Receptor beta antibody (orb575669)

Description	Rabbit polyclonal antibody to Estrogen Receptor beta
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human ESR2
Target	ESR2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SYVDSHHEYPMATFYSPAVMNYSIPSNVTNLEGGPGRQTTPNVLWPTPG
Clonality	Polyclonal
MW	56kDa
Uniprot ID	Q92731
NCBI	NP_001428
Expiration Date	12 months from date of receipt.