


Product Datasheet

EBF1 antibody (orb577182)

Description	Rabbit polyclonal antibody to EBF1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards (50ug) human EBF1.
Target	EBF1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	LPSNCSSSSGIFSFS PANMVSAVKQKSAFAPVVRPQTSPPTCTSTNGNS
Clonality	Polyclonal
Purity	Affinity Purified
MW	64 kDa
Uniprot ID	Q9UH73
Expiration Date	12 months from date of receipt.



WB Suggested Anti-EBF1 Antibody Titration...



Immunohistochemistry with Thymus tissue ...