

Product Datasheet

Caveolin-1 (phospho-Tyr14) antibody (orb5773)

Description	Caveolin-1 (phospho-Tyr14) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Caveolin-1 around the phosphorylation site of Tyr14 (HL(p-Y)TV)
Target	Phospho-Caveolin-1 (Tyr14)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	Q03135
Expiration Date	12 months from date of receipt.