

Product Datasheet

CACH2 antibody (orb5775)

Description	CACH2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CACH2 (351-450/2221 aa)
Target	CACNA1C
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	239 kDa
Uniprot ID	Q13936
Expiration Date	12 months from date of receipt.