


Product Datasheet


Doublecortin antibody (orb578005)

Description	Rabbit polyclonal antibody to Doublecortin
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human DCX
Target	DCX
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PEKFRYAQDDFSLDENECRVMKGNPSATAGPKASPTPQKTSKSPGPMRR
Clonality	Polyclonal
MW	40kDa
Uniprot ID	O43602
NCBI	NP_835365
Expiration Date	12 months from date of receipt.



See PC 1 Data and customer Feedback for more information

Immunohistochemical
staining of Mouse
sp...



Western blot analysis
of HepG2 cell lysa...