

Product Datasheet

CBL2 (phospho-Tyr731) antibody (orb5781)

Description	CBL2 (phospho-Tyr731) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CBL2 around the phosphorylation site of Tyr731 (CT(p-Y)EA)
Target	phospho-CBL (Tyr731)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100 kDa
Uniprot ID	P22681
Expiration Date	12 months from date of receipt.