

Product Datasheet

ALAS2 antibody (orb578150)

Description	Rabbit polyclonal antibody to ALAS2
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Goat, Guinea pig, Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ALAS2
Target	ALAS2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P22557
Expiration Date	12 months from date of receipt.