

## Product Datasheet

### PLEKHA4 antibody (orb578525)

|                            |                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PLEKHA4                                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                                                   |
| <b>Reactivity</b>          | Animal, Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rat                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                             |
| <b>Tested Applications</b> | WB                                                                                                                                       |
| <b>Immunogen</b>           | The immunogen is a synthetic peptide directed towards the middle region of human PLEKHA4                                                 |
| <b>Target</b>              | PLEKHA4                                                                                                                                  |
| <b>Preservatives</b>       | Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                                |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                   |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                                                |
| <b>Storage</b>             | For short term use, store at 2-8C up to 1 week.<br>For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                                                               |
| <b>Uniprot ID</b>          | <b>Q9H4M7</b>                                                                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                          |