


Product Datasheet

CD19 (Phospho-Tyr531) antibody (orb5789)

Description	Rabbit polyclonal to CD19 (Phospho-Tyr531). CD19 a cell surface molecule which assembles with the an...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...