



Product Datasheet

SLC10A5 antibody (orb578928)

Description	Rabbit polyclonal antibody to SLC10A5
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human SLC10A5
Target	SLC10A5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	Batch dependent within range: 100 ul at 0.5 - 1 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GYSFAKVCTLPVCKTVAIESGMLNSFLALAVIQLSFPQSKANLASVAP
Clonality	Polyclonal
MW	48kDa
Uniprot ID	Q5PT55
NCBI	NP_001010893
Expiration Date	12 months from date of receipt.



Western blot analysis of human Jurkat cells...



Immunohistochemical staining of human lung...