



Product Datasheet

ADCY7 antibody (orb5795)

Description	ADCY7 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human ADCY7 (901-1080/1080 aa)
Target	ADCY7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120 kDa
Uniprot ID	P51828
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Hela cells
U...



Immunohistochemical staining of rat colo...