



Product Datasheet

ADCY3 antibody (orb5798)

Description	ADCY3 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ADCY3 (301-400/1144 aa)
Target	ADCY3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	126 kDa
Uniprot ID	O60266
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Immunohistochemical staining of rat brai...