



Product Datasheet

CD44v4 antibody (orb5801)


Description	CD44v4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human CD44v4 (141-240/361 aa)
Target	CD44v4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P16070
Dilution Range	ICC=1:100, Flow-Cyt=1µg /Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...



Western blot analysis of extracts from b...



Flow cytometric analysis of Mouse spleen...