


Product Datasheet

CD158b2 antibody (orb5811)

Description	CD158b2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human CD158b2 (42-140/341 aa)
Target	KIR2DL3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22/38 kDa
Uniprot ID	P43628
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse CO...