

Product Datasheet

LRRFIP2 antibody (orb581145)

Description	Rabbit polyclonal antibody to LRRFIP2
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rat, Yeast, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	LRRFIP2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KQYAENYTRPSSRNSASATTPLSGNSSRRGSGDTSSLIDPDTSLSELRES
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q9Y608
NCBI	NP_060194
Expiration Date	12 months from date of receipt.