



Product Datasheet

NUP50 antibody (orb581651)


Description	Rabbit polyclonal antibody to NUP50
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human NUP50
Target	NUP50
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	100 µl at 0.5 - 1 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TTQSKPVSSPFPTKPLEGQAEGDSGECKGGDEEENDEPPKVVVTEVKEED
Clonality	Polyclonal
MW	50kDa
Uniprot ID	Q9UKX7
NCBI	NP_009103
Expiration Date	12 months from date of receipt.



Western blot analysis of human hela cell...



WB analysis of HepG2 tissue using NUP50 ...



Immunofluorescence analysis of human hel...