

Product Datasheet

Nudt16 antibody (orb581885)

Description	Rabbit polyclonal antibody to Nudt16
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Target	Nudt16
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VLMQMRFDGRLGFPGGFVDAQDSCLEDGLNRELREELGEAMSAFRVERSD
Clonality	Polyclonal
MW	23kDa
Uniprot ID	Q6P3D0
NCBI	NP_083661
Expiration Date	12 months from date of receipt.