


Product Datasheet

cdc25C (phospho-Thr48) antibody (orb5826)

Description	cdc25C (phospho-Thr48) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human cdc25C around the phosphorylation site of Thr48 (PR(p-T)PV)
Target	Phospho-cdc25C (Thr48)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	P30307
Dilution Range	IHC-P=1:100-500, ICC=1:100, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of Hela cell using Phospho-...