

Product Datasheet

C3orf19 antibody (orb583367)

Description	Rabbit polyclonal antibody to C3orf19
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Guinea pig, Human, Mouse, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human C3orf19
Target	CCDC174
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RQQWEEEEEREALKRPMGVPVHYEDIRENEARQLGVGYFAFARDKELRNKQM
Clonality	Polyclonal
MW	54kDa
Uniprot ID	Q6PII3
NCBI	NP_057558
Expiration Date	12 months from date of receipt.