

Product Datasheet

Gpbp1l1 antibody (orb583611)

Description	Rabbit polyclonal antibody to Gpbp1l1
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen for Anti-Gpbp1l1 antibody is: synthetic peptide directed towards the C-terminal of Mouse GPBL1
Target	Gpbp1l1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AGAGLHSPKESPSSTTPPIEISSRLTKLTRRTDRKSEFLKTLKDERNG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	Q6NZP2
NCBI	XP_006503548.1
Expiration Date	12 months from date of receipt.