


Product Datasheet

CEBP alpha (phospho-Ser21) antibody (orb5840)

Description	CEBP alpha (phospho-Ser21) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CDBP alpha around the phosphorylation site of Ser21 (LQ(p-S)PP)
Target	Phospho-CEBP alpha (Ser21)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	P49715
Dilution Range	ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



ICC analysis of HepG2 cell using Phospho...