

Product Datasheet

ZCCHC7 antibody (orb584714)

Description	Rabbit polyclonal antibody to ZCCHC7
Species/Host	Rabbit
Reactivity	Animal, Human, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human ZCCHC7
Target	ZCCHC7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KNRNWEKHKADRHREVDDEFPRGPKTYSSPGSFKTQKPSKPFHRSSHYH
Clonality	Polyclonal
MW	63kDa
Uniprot ID	Q8N3Z6
NCBI	NP_115602
Expiration Date	12 months from date of receipt.