


Product Datasheet

CTPS2 antibody (orb585338)

Description	Rabbit polyclonal antibody to CTPS2
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rabbit, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Target	CTPS2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AKRENFCAIHVSLVPQLSATGEQTKPTQNSVRALRGLGLSPDLIVCRSS
Clonality	Polyclonal
MW	64kDa
Uniprot ID	Q9NRF8
NCBI	NP_787055
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human He...



Western blot analysis of HeLa Whole Cell...