


Product Datasheet

c-FOS(phospho-Thr232) antibody (orb5856)

Description	c-FOS(phospho-Thr232) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse c-Fos around the phosphorylation site of Thr232 (AS(p-T)PE)
Target	Phospho-c-Fos (Thr232)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P01100
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western
Blot using c-
Fos
(phospho-
Thr232...