



Product Datasheet

Agpat4 antibody (orb5860)


Description	Agpat4 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Agpat4 (231-330/378 aa)
Target	AGPAT4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q9NRZ5
Expiration Date	12 months from date of receipt.



IHC-P of mouse carcinoma tissue using Ag...



IHC-P of rat pancreatic tissue(Agpat4 an...



Western blot analysis of rat pancreatic ...