

Product Datasheet

c-Jun (phospho-Ser73) antibody [Out of stock] (orb5865)

Description	c-Jun (phospho-Ser73) antibody [Out of stock]
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Jun around the phosphorylation site of Ser73
Target	phospho-c-Jun (Ser73)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05412
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-500, IHC-F: 1:100-500, IF: 1:100-500, Flow-Cyt: 1µg /Test, ELISA: 1:5000-10000
Expiration Date	12 months from date of receipt.