

Product Datasheet

HSPA12A antibody (orb586612)

Description	Rabbit polyclonal antibody to HSPA12A
Species/Host	Rabbit
Reactivity	Animal, Bovine, Canine, Goat, Guinea pig, Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human HSPA12A
Target	HSPA12A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KALEIFAYALQYFKEQALKELSDQAGSEFENS DVRWVITVPAIWKQPAKQ
Clonality	Polyclonal
MW	74kDa
Uniprot ID	O43301
NCBI	NP_079291
Expiration Date	12 months from date of receipt.