


Product Datasheet

PCNA (phospho-Tyr211) antibody (orb5867)

Description	PCNA (phospho-Tyr211) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PCNA around the phosphorylation site of Tyr211 (LR(p-Y)LN)
Target	phosphos-PCNA (Tyr211)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29 kDa
Uniprot ID	P12004
Dilution Range	ICC=1:100, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...