

Product Datasheet

TIGD5 antibody (orb587168)

Description	Rabbit polyclonal antibody to TIGD5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human TIGD5
Target	TIGD5
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SVC R D F G V P G G T L R G W L K D E P K L R W F L E Q L G G E V G T Q R K K M R L A N E E E I D
Clonality	Polyclonal
MW	70 kDa
Uniprot ID	Q53EQ6
NCBI	NP_116251.4
Expiration Date	12 months from date of receipt.