

Product Datasheet

RAD51D antibody (orb588429)

Description	Rabbit polyclonal antibody to RAD51D
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen for Anti-RAD51D antibody is: synthetic peptide directed towards the Middle region of Human RA51D
Target	RAD51D
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SPLGGGQREGALMMQLARELKTLDLGMVVVTNHITRDRDSGRLKP
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	O75771
NCBI	NP_598332.1
Expiration Date	12 months from date of receipt.