

Product Datasheet

RPGRIP1L antibody (orb589840)

Description	Rabbit polyclonal antibody to RPGRIP1L
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle terminal region of human RPGRIP1L
Target	RPGRIP1L
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VLAPRPKPRQRLTPVDKKVSFVDIMPHQSDETSPPPEDRKEISPEVEHIP
Clonality	Polyclonal
MW	144 kDa
Uniprot ID	Q68CZ1
NCBI	NP_001121369.1
Expiration Date	12 months from date of receipt.